

Product Datasheet

Ac-Hirudin (55-65) (Desulfated) (orb364145)

Description	Peptide to Ac-Hirudin (55-65) (Desulfated)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ac-Asp-Phe-Glu-Glu-Ile-Pro-Glu-Glu-Tyr-Leu-Gln
Formula	$C_{66}H_{92}N_{12}O_{25}$
Purity	> 95%
Source	Synthetic
MW	1453.53
CAS Number	113274-57-0
Expiration Date	12 months from date of receipt.

