

Product Datasheet

Papain Substrate Colorimetric (orb364153)

Description	Peptide to Papain Substrate Colorimetric
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ac-Phe-Gly-pNA
Formula	C ₁₉ H ₂₀ N ₄ O ₅
Purity	> 95%
Source	Synthetic
MW	384.40
CAS Number	34336-99-7
Expiration Date	12 months from date of receipt.

