

Product Datasheet

Adipokinetic Hormone II (*Locusta migratoria*) (orb364170)

Description	Peptide to Adipokinetic Hormone II (<i>Locusta migratoria</i>)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Pyr-Leu-Asn-Phe-Ser-Ala-Gly-Trp-NH ₂
Formula	C ₄₃ H ₅₇ N ₁₁ O ₁₁
Purity	> 95%
Source	Synthetic
MW	904.02
Expiration Date	12 months from date of receipt.