

Product Datasheet

Alarelin Acetate (orb364187)

Description	Peptide to Alarelin Acetate
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	pGlu-His-Trp-Ser-Tyr-D-Ala-Leu-Arg-Pro-NHC2H5
Formula	$C_{56}H_{78}N_{16}O_{12}$
Purity	> 95%
Source	Synthetic
MW	1167.3
CAS Number	52435-06-0
Expiration Date	12 months from date of receipt.