

Product Datasheet

Allatotropin (orb364194)

Description	Peptide to Allatotropin
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Gly-Phe-Lys-Asn-Val-Glu-Met-Met-Thr-Ala-Arg-Gly-Phe-NH ₂ (Disulfide bridge Cys2-Cys7)
Formula	C ₆₅ H ₁₀₃ N ₁₉ O ₁₇ S ₂
Purity	> 95%
Source	Synthetic
MW	1486.80
CAS Number	120928-88-3
Expiration Date	12 months from date of receipt.

