

Product Datasheet

Alloferon 1 (orb364195)

Description	Peptide to Alloferon 1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	His-Gly-Val-Ser-Gly-His-Gly-Gln-His-Gly-Val-His-Gly
Formula	$C_{52}H_{76}N_{22}O_{16}$
Purity	> 95%
Source	Synthetic
MW	1265.30
CAS Number	347884-61-1
Expiration Date	12 months from date of receipt.

