

Product Datasheet

Amyloid Bri Protein (1-34) (orb364229)

Description	Peptide to Amyloid Bri Protein (1-34)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Pyr-Ala-Ser-Asn-Cys-Phe-Ala-Ile-Arg-His-Phe-Glu-Asn-Lys-Phe-Ala-Val-Glu-Thr-Leu-Ile-Cys-Ser-Arg-Thr-Val-Lys-Lys-Asn-Ile-Ile-Glu-Glu-Asn (Disulfide bridge Cys5-Cys22)
Formula	$C_{173}H_{273}N_{49}O_{52}S_2$
Purity	> 95%
Source	Synthetic
MW	3935.55
Expiration Date	12 months from date of receipt.