

Product Datasheet

Angiotensin (1-5) (orb364242)

Description	Peptide to Angiotensin (1-5)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Asp-Arg-Val-Tyr-Ile
Formula	$C_{30}H_{48}N_8O_9$
Purity	> 95%
Source	Synthetic
MW	664.77
CAS Number	58442-64-1
Expiration Date	12 months from date of receipt.

