


Product Datasheet

Angiotensin I-Converting Enzyme Substrate (orb364247)

Description	Peptide to Angiotensin I-Converting Enzyme Substrate
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Bz-Phe-Ala-Pro
Formula	$C_{24}H_{27}N_3O_5$
Purity	> 95%
Source	Synthetic
MW	437.49
CAS Number	69677-91-4
Expiration Date	12 months from date of receipt.



Angiotensin
I-Converting
Enzyme
Substrat...