

Product Datasheet

Angiotensinogen peptide (orb364251)

Description	Angiotensinogen peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Asp-Arg-Val-Tyr-Ile-His-Pro-Phe-His-Leu-Val-Ile-His
Formula	$C_{79}H_{116}N_{22}O_{17}$
Purity	> 95%
Source	Synthetic
MW	1645.90
CAS Number	82048-97-3
Expiration Date	12 months from date of receipt.

