

Product Datasheet

Angiotensinogen peptides (orb364252)

Description	Angiotensinogen peptides
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Asp-Arg-Val-Tyr-Ile-His-Pro-Phe-His-Leu-Leu-Val-Tyr-Ser
Formula	$C_{85}H_{123}N_{21}O_{20}$
Purity	> 95%
Source	Synthetic
MW	1759.10
CAS Number	64315-16-8
Expiration Date	12 months from date of receipt.

