

Product Datasheet

Antifreeze Polypeptide (orb364261)

Description	Peptide to Antifreeze Polypeptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Asp-Thr-Ala-Ser-Asp-Ala-Ala-Ala-Ala-Ala-Ala-Leu-Thr-Ala-Ala-Asn-Ala-Lys-Ala-Ala-Ala-Glu-Leu-Thr-Ala-Ala-Asn-Ala-Ala-Ala-Ala-Ala-Ala-Thr-Ala-Arg
Formula	$C_{133}H_{225}N_{43}O_{51}$
Purity	> 95%
Source	Synthetic
MW	3242.53
Expiration Date	12 months from date of receipt.