

Product Datasheet

Anti-Inflammatory Peptide 3 (orb364264)

Description	Peptide to Anti-Inflammatory Peptide 3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Met-Gln-Met-Asn-Lys-Val-Leu-Asp-Ser
Formula	$C_{43}H_{76}N_{12}O_{15}S_2$
Purity	> 95%
Source	Synthetic
MW	1065.28
CAS Number	118850-73-0
Expiration Date	12 months from date of receipt.

