

Product Datasheet

Apelin-15 (63-77) (orb364272)

Description	Peptide to Apelin-15 (63-77)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Arg-Arg-Gln-Arg-Pro-Arg-Leu-Ser-His-Lys- Gly-Pro-Met-Pro-Phe
Formula	$C_{81}H_{135}N_{31}O_{18}S$
Purity	> 95%
Source	Synthetic
MW	1863.24
Expiration Date	12 months from date of receipt.