

Product Datasheet

Aquaporin 2 Peptides (255-271) (orb364276)

Description	Aquaporin 2 Peptides (255-271)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Gln-Ser-Val-Glu-Leu-His-Ser-Pro-Gln-Ser-Leu-Pro-Arg-Gly-Thr-Lys-Ala
Formula	$C_{78}H_{131}N_{25}O_{26}$
Purity	> 95%
Source	Synthetic
MW	1835.07
Expiration Date	12 months from date of receipt.