

Product Datasheet

ASIF (1-35) (orb364280)

Description	Peptide to ASIF (1-35) (Bovine)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ala-Leu-Arg-Gly-Pro-Lys-Met-Met-Arg-Asp-Ser-Gly-Cys-Phe-Gly-Arg-Arg-Leu-Asp-Arg-Ile-Gly-Ser-Leu-Ser-Gly-Leu-Gly-Cys-Asn-Val-Leu-Arg-Arg-Tyr (Disulfide bridge Cys13-Cys29)
Formula	$C_{164}H_{278}N_{58}O_{45}S_4$
Purity	> 95%
Source	Synthetic
MW	3910.64
CAS Number	120249-06-1
Expiration Date	12 months from date of receipt.