

Product Datasheet

Angiotensinogen receptor peptide (orb364281)

Description	Angiotensinogen receptor peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ala-Tyr-Glu-Ile-Gln-Lys-Asn-Lys-Pro-Arg-Asn-Asp-Asp
Formula	$C_{67}H_{107}N_{21}O_{24}$
Purity	> 95%
Source	Synthetic
MW	1590.73
Expiration Date	12 months from date of receipt.