

Product Datasheet

Band 3 Protein (547-553) (orb364292)

Description	Peptide to Band 3 Protein (547-553) (Human)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	His-Pro-Leu-Gln-Lys-Thr-Tyr
Formula	$C_{41}H_{63}N_{11}O_{11}$
Purity	> 95%
Source	Synthetic
MW	886.02
CAS Number	167095-71-8
Expiration Date	12 months from date of receipt.

