

Product Datasheet

Band 3 Protein (824-829) (orb364293)

Description	Peptide to Band 3 Protein (824-829) (Human)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Tyr-Val-Lys-Arg-Val-Lys
Formula	$C_{37}H_{65}N_{11}O_8$
Purity	> 95%
Source	Synthetic
MW	791.99
CAS Number	158475-15-1
Expiration Date	12 months from date of receipt.

