

Product Datasheet

BAX peptide (orb364294)

Description	BAX peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Lys-Lys-Leu-Ser-Glu-Cys-Leu-Lys-Arg-Ile-Gly-Asp-Glu-Leu-Asp-Ser
Formula	$C_{77}H_{136}N_{22}O_{27}S$
Purity	> 95%
Source	Synthetic
MW	1834.13
Expiration Date	12 months from date of receipt.