

Product Datasheet

Beta II probe (orb364298)

Description	Peptide to Beta II probe
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ser-Ser-Leu-Asp-Leu-Ser-Gln-Phe-Pro-Met-Thr-Ala-Ser-Phe-Leu-Arg-Glu-Ser-Arg
Formula	$C_{94}H_{150}N_{26}O_{31}S$
Purity	> 95%
Source	Synthetic
MW	2172.46
Expiration Date	12 months from date of receipt.

