

Product Datasheet

Beta-Amyloid (10-35) (orb364301)

Description	Peptide to Beta-Amyloid (10-35)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Tyr-Glu-Val-His-His-Gln-Lys-Leu-Val-Phe-Phe-Ala-Glu-Asp-Val-Gly-Ser-Asn-Lys-Gly-Ala-Ile-Ile-Gly-Leu-Met
Formula	$C_{133}H_{204}N_{34}O_{37}S$
Purity	> 95%
Source	Synthetic
MW	2903.38
Expiration Date	12 months from date of receipt.