

Product Datasheet

Beta-Amyloid (2-40) (orb364329)

Description	Peptide to Beta-Amyloid (2-40)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ala-Glu-Phe-Arg-His-Asp-Ser-Gly-Tyr-Glu-Val-His-His-Gln-Lys-Leu-Val-Phe-Phe-Ala-Glu-Asp-Val-Gly-Ser-Asn-Lys-Gly-Ala-Ile-Ile-Gly-Leu-Met-Val-Gly-Gly-Val-Val
Formula	$C_{190}H_{290}N_{52}O_{55}S$
Purity	> 95%
Source	Synthetic
MW	4214.81
Expiration Date	12 months from date of receipt.

