

Product Datasheet

Beta-Amyloid (96-110) analog (orb364345)

Description	Peptide to Beta-Amyloid (96-110) analog
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ac-Asn-Trp-Cys-Lys-Arg-Gly-Arg-Lys-Gln-Cys-Lys-Thr-His-Pro-His-NH ₂ (Disulfide bridges Cys98-Cys105)
Formula	C ₈₁ H ₁₂₈ N ₃₂ O ₁₉ S ₂
Purity	> 95%
Source	Synthetic
MW	1918.25
Expiration Date	12 months from date of receipt.

