

Product Datasheet

scn antibody (FITC) (orb36435)

Description	scn antibody (FITC)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	FITC
Immunogen	Recombinant Staphylococcus aureus Staphylococcal complement inhibitor protein (32-116AA)
Target	scn
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q99SU9
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.