

Product Datasheet

Beta-Lipotropin (88-91) (orb364366)

Description	Peptide to Beta-Lipotropin (88-91)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Lys-Lys-Gly-Glu
Formula	$C_{19}H_{36}N_6O_7$
Purity	> 95%
Source	Synthetic
MW	460.53
CAS Number	72189-84-5
Expiration Date	12 months from date of receipt.

