

Product Datasheet

bFGF Inhibitory Peptide (orb364372)

Description	Peptide to bFGF Inhibitory Peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ala-Pro-Ser-Gly-His-Tyr-Lys-Gly
Formula	$C_{36}H_{53}N_{11}O_{11}$
Purity	> 95%
Source	Synthetic
MW	815.89
CAS Number	152051-60-0
Expiration Date	12 months from date of receipt.

