

Product Datasheet

Bid BH3 (orb364374)

Description	Peptide to Bid BH3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Arg-Asn-Ile-Ala-Arg-His-Leu-Ala-Gln-Val-Gly-Asp-Ser-Met-Asp-Arg
Formula	$C_{74}H_{127}N_{29}O_{24}S$
Purity	> 95%
Source	Synthetic
MW	1839.08
Expiration Date	12 months from date of receipt.