

Product Datasheet

Bid BH3 Peptide (orb364375)

Description	Peptide to Bid BH3 Peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Glu-Asp-Ile-Ile-Arg-Asn-Ile-Ala-Arg-His-Leu-Ala-Gln-Val-Gly-Asp-Ser-Met-Asp-Arg
Formula	$C_{95}H_{161}N_{33}O_{32}S$
Purity	> 95%
Source	Synthetic
MW	2309.61
Expiration Date	12 months from date of receipt.