

Product Datasheet

Bid BH3-r8 (orb364376)

Description	Peptide to Bid BH3-r8
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	D-Arg-D-Arg-D-Arg-D-Arg-D-Arg-D-Arg-D-Arg-Gly-Glu-Asp-Ile-Ile-Arg-Asn-Ile-Ala-Arg-His-Leu-Ala-Gln-Val-Gly-Asp-Ser-Met-Asp-Arg
Formula	$C_{145}H_{260}N_{66}O_{41}S$
Purity	> 95%
Source	Synthetic
MW	3616.17
Expiration Date	12 months from date of receipt.