

Product Datasheet

Bid BH3-r9 (orb364377)

Description	Peptide to Bid BH3-r9
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Arg-Arg-Arg-Arg-Arg-Arg-Arg-Arg-Gly-Glu-Asp-Ile-Ile-Arg-Asn-Ile-Ala-Arg-His-Leu-Ala-Gln-Val-Gly-Asp-Ser-Met-Asp-Arg
Formula	$C_{151}H_{272}N_{70}O_{42}S$
Purity	> 95%
Source	Synthetic
MW	3772.36
Expiration Date	12 months from date of receipt.