

Product Datasheet

Vasopressin peptide (Biotin) (orb364383)

Description	Peptide to Biotin-(Arg8)-Vasopressin
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Biotin-Cys-Tyr-Phe-Gln-Asn-Cys-Pro-Arg-Gly-NH ₂ (Disulfide bridge Cys1-Cys6)
Formula	C ₅₆ H ₇₉ N ₁₇ O ₁₄ S ₂
Purity	> 95%
Source	Synthetic
MW	1310.55
Expiration Date	12 months from date of receipt.