

Product Datasheet

Biotin-[Glu1]-Gastrin Receptor I (phosphorylated) [Out of stock] (orb364385)

Description	Biotin-[Glu1]-Gastrin Receptor I (phosphorylated)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Biotin-Glu-Gly-Pro-Trp-Leu-Glu-Glu-Glu-Glu-Glu-Ala-pTyr-Gly-Trp-Met-Asp-Phe-NH ₂
Formula	C ₁₀₇ H ₁₄₁ N ₂₂ O ₃₇ PS ₂
Purity	> 95%
Source	Synthetic
MW	2422.53
Expiration Date	12 months from date of receipt.

