

Product Datasheet

Biotin-Atrial Natriuretic Peptide (1-28) (orb364392)

Description	Peptide to Biotin-Atrial Natriuretic Peptide (1-28) (Human, Porcine)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Biotin-Ser-Leu-Arg-Arg-Ser-Ser-Cys-Phe-Gly-Gly-Arg-Met-Asp-Arg-Ile-Gly-Ala-Gln-Ser-Gly-Leu-Gly-Cys-Asn-Ser-Phe-Arg-Tyr (Disulfide bridge Cys7-Cys23)
Formula	$C_{137}H_{217}N_{47}O_{41}S_4$
Purity	> 95%
Source	Synthetic
MW	3308.80
Expiration Date	12 months from date of receipt.