

Product Datasheet

BNP-26 (orb364431)

Description	Peptide to BNP-26 (Porcine)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Asp-Ser-Gly-Cys-Phe-Gly-Arg-Arg-Leu-Asp-Arg-Ile-Gly-Ser- Leu-Ser-Gly-Leu-Gly-Cys-Asn-Val-Leu-Arg-Arg-Tyr (Disulfide bridge Cys4-Cys5)
Formula	$C_{120}H_{198}N_{42}O_{36}S_2$
Purity	> 95%
Source	Synthetic
MW	2869.30
CAS Number	114547-28-3
Expiration Date	12 months from date of receipt.

