

Product Datasheet

BPDEtide (orb364439)

Description	Peptide to BPDEtide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Arg-Lys-Ile-Ser-Ala-Ser-Glu-Phe-Asp-Arg-Pro-Leu-Arg
Formula	$C_{68}H_{115}N_{23}O_{20}$
Purity	> 95%
Source	Synthetic
MW	1574.82
Expiration Date	12 months from date of receipt.

