

Product Datasheet

cAMP Dependent PK Inhibitor (5-24) (orb364468)

Description	Peptide to cAMP Dependent PK Inhibitor (5-24)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Thr-Thr-Tyr-Ala-Asp-Phe-Ile-Ala-Ser-Gly-Arg-Thr-Gly-Arg-Arg-Asn-Ala-Ile-His-Asp
Formula	$C_{94}H_{108}N_{32}O_{31}$
Purity	> 95%
Source	Synthetic
MW	2222.40
Expiration Date	12 months from date of receipt.