

Product Datasheet

cAMP Dependent Protein Kinase Substrate (orb364469)

Description	Peptide to cAMP Dependent Protein Kinase Substrate
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Arg-Arg-Lys-Ala-Ser-Gly-Pro
Formula	$C_{31}H_{58}N_{14}O_9$
Purity	> 95%
Source	Synthetic
MW	770.90
CAS Number	65189-70-0
Expiration Date	12 months from date of receipt.

