

Product Datasheet

CAP (orb364470)

Description	Peptide to CAP (Rat)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Asp-Met-Ala-Lys-Asp-Leu-Glu-Thr-Asn-His-His-Pro-Tyr-Phe-Gly-Asn
Formula	$C_{82}H_{117}N_{23}O_{27}S$
Purity	> 95%
Source	Synthetic
MW	1889.05
Expiration Date	12 months from date of receipt.

