

Product Datasheet

Chorionic Gonadotropin-Beta (109-145) (orb364514)

Description	Peptide to Chorionic Gonadotropin-Beta (109-145) (Human)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Thr-Cys-Asp-Asp-Pro-Arg-Phe-Gln-Asp-Ser-Ser-Ser-Ser-Lys-Ala-Pro-Pro-Pro-Ser-Leu-Pro-Ser-Pro-Ser-Arg-Leu-Pro-Gly-Pro-Ser-Asp-Thr-Pro-Ile-Leu-Pro-Gln
Formula	$C_{167}H_{264}N_{46}O_{58}S_1$
Purity	> 95%
Source	Synthetic
MW	3876.22
CAS Number	72979-70-5
Expiration Date	12 months from date of receipt.