

Product Datasheet

C-Reactive Protein (77-82) (orb364531)

Description	Peptide to C-Reactive Protein (77-82)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Val-Gly-Gly-Ser-Glu-Ile
Formula	$C_{23}H_{40}N_6O_{10}$
Purity	> 95%
Source	Synthetic
MW	560.61
Expiration Date	12 months from date of receipt.