

Product Datasheet

CSH103 (orb364536)

Description	Peptide to CSH103
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ala-Asp-Ala-Gln-His-Ala-Thr-Pro-Pro-Lys-Lys-Lys-Arg-Lys-Val-Glu-Asp-Pro-Lys-Asp-Phe
Formula	$C_{106}H_{172}N_{32}O_{32}$
Purity	> 95%
Source	Synthetic
MW	2406.80
CAS Number	135546-44-0
Expiration Date	12 months from date of receipt.

