

Product Datasheet

Dendroaspis Natriuretic Peptide (orb364554)

Description	Peptide to Dendroaspis Natriuretic Peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Glu-Val-Lys-Tyr-Asp-Pro-Cys-Phe-Gly-His-Lys-Ile-Asp-Arg-Ile-Asn-His-Val-Ser-Asn-Leu-Gly-Cys-Pro-Ser-Leu-Arg-Asp-Pro-Arg-Pro-Asn-Ala-Pro-Ser-Thr-Ser-Ala (Disulfide bridge Cys7-Cys23)
Formula	$C_{180}H_{282}N_{56}O_{56}S_2$
Purity	> 95%
Source	Synthetic
MW	4190.72
CAS Number	255721-52-9
Expiration Date	12 months from date of receipt.