

Product Datasheet

Doc-6 (130-145) (orb364557)

Description	Peptide to Doc-6 (130-145)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Gly-Val-Gln-Arg-Glu-Gln-Asn-Glu-Arg-Phe-Asn-Val-Tyr-Leu-Met-Pro
Formula	$C_{86}H_{134}N_{26}O_{26}S$
Purity	> 95%
Source	Synthetic
MW	1980.25
Expiration Date	12 months from date of receipt.