

Product Datasheet

Dok-1 (348-359) (orb364558)

Description	Peptide to Dok-1 (348-359)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Leu-Leu-Lys-Ala-Lys-Leu-Thr-Asp-Pro-Lys-Glu-Asp
Formula	$C_{61}H_{107}N_{15}O_{20}$
Purity	> 95%
Source	Synthetic
MW	1370.62
Expiration Date	12 months from date of receipt.