

Product Datasheet

Dok-4 (130-145) (orb364559)

Description	Peptide to Dok-4 (130-145)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Gly-Val-Gln-Cys-Glu-Gln-Thr-Asp-Arg-Phe-Asn-Val-Phe-Leu-Leu-Pro
Formula	$C_{83}H_{128}N_{22}O_{25}S$
Purity	> 95%
Source	Synthetic
MW	1866.10
Expiration Date	12 months from date of receipt.