

Product Datasheet

Dok-4 (263-275) (orb364560)

Description	Peptide to Dok-4 (263-275)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Leu-Pro-Arg-Ser-Ala-Tyr-Trp-His-His-Ile-Thr-Gly-Ser
Formula	$C_{70}H_{101}N_{21}O_{18}$
Purity	> 95%
Source	Synthetic
MW	1524.72
Expiration Date	12 months from date of receipt.