

Product Datasheet

Dynorphin A (1-6) (orb364571)

Description	Peptide to Dynorphin A (1-6) (Porcine)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Tyr-Gly-Gly-Phe-Leu-Arg
Formula	$C_{34}H_{49}N_9O_8$
Purity	> 95%
Source	Synthetic
MW	711.83
CAS Number	75106-70-6
Expiration Date	12 months from date of receipt.

