

Product Datasheet

Dynorphin A (1-8) (orb364573)

Description	Peptide to Dynorphin A (1-8) (Porcine)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Tyr-Gly-Gly-Phe-Leu-Arg-Arg-Ile
Formula	$C_{46}H_{72}N_{14}O_{10}$
Purity	> 95%
Source	Synthetic
MW	981.17
CAS Number	75790-53-3
Expiration Date	12 months from date of receipt.

