

Product Datasheet

Dynorphin A (2-12) (orb364575)

Description	Peptide to Dynorphin A (2-12) (Porcine)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Gly-Gly-Phe-Leu-Arg-Arg-Ile-Arg-Pro-Lys-Leu
Formula	$C_{60}H_{105}N_{21}O_{12}$
Purity	> 95%
Source	Synthetic
MW	1312.64
Expiration Date	12 months from date of receipt.

